

POSY2 Quick Reference

Instrument computer: **posy2.pns.anl.gov** (VAX/VMS)

Username: POSY

Password: *see a member of the POSY2 staff*

Common commands:

newrun set up a new run

start or **start 12540-12542** or **start 12540,12541,12542**

stoprun stop the run and save the data

quitrn stop the run without saving the data

run align run program that aligns the sample

ctl 12540 change control parameters of run 12540

prun control 12540 check control parameters of run 12540

schedule 2(12540)/108000P schedule run 12540 to run for 2 cycles (= 2hours) of 108000 pulses

Looking at data:

typ lh look at data actively being collected (live histogram)

typ fh choose to look at data of a previous run (filed histogram)

runn 12540 specify run 12540 as the data you want to look at

hst x look at data as a function of position channel **x** on the detector

tpt 1,256 select to look at the sum of time channels 1,256

dis display plot

cvs close plot

hst t look at data as a function of time channel **t** (time = wavelength)
 $\text{lambda (\AA)} = 5.10\text{e-}4 * t (\text{microsec})$

xpt 20,30 select to look at the sum of position channels 20-30

rdp read number of measured pulses

integrate 1,10 integrate over channels 1-10 as shown in the plot

[@clearmem](#) clear the computer memory buffer while a measurement is running (the data will be lost)

[ymax](#) or [ymin](#) set the min or max value for the y axis.

[ysc](#) change the yscale to the newly defined ymax or ymin, or revert back to the auto scaling

[dx 10,230](#) set the x axis min and max to 10 and 230, respectively

[prun/live status](#) check status of live run

[prun status 12540](#) check final status of run 12540 (info on # of measured pulses in file, title etc.)

[prun header 12540](#) Get header info out of runfile of run 12540, such as # of measured pulses in file, title, total number of monitor counts, monitor counts per channel, total number of detector counts

[dintall](#) if you get the message that a display can't be opened because the "graphics are busy" type this command

Ancillary equipment:

[setang](#) change reflection angle

[translate](#) change translation of sample perpendicular to beam

[angcal](#) check the angle calibration

[rdt](#) read temperature

[temtur](#) change temperature

Logging in as the instrument account on the LINUX system at the instrument

Username: *posy2* *case sensitive !!!*

Password: *see a member of the POSY2 staff* *case sensitive !!!*

Data analysis: **aragorn.pns.anl.gov**, **boromir.pns.anl.gov** or **gimli.pns.anl.gov** (LINUX)

Username: *your IPNS LINUX Cluster username* *case sensitive !!!*

Password: *your IPNS LINUX Cluster password* *case sensitive !!!*

Run IDL program:

[posyidl2](#)

Select [POSY2](#) button

Select menu [Data Analysis](#)

[Analyze Reflectivity Data](#)

Select menu [Data Analysis](#)

Merge Reflectivity Files (automatic merging)
Or Manually Merge Reflectivity Files (manual merging)
Select menu Plotting
Analyzed Data
Select menu Calculation
Layers
Or Refit
Also:
Select menu Plotting
Contour (make 2D contour plot of angle vs wavelength or qx vs qz)

Remote login to VAX/VMS computers via LINUX systems:

1) Log into **dasmmain.pns.anl.gov** using same username and password as for **gimli.pns.anl.gov**
type: rlogin -l posy posy2 to log into instrument computer **posy2.pns.anl.gov**